OXYCYCLE 2

Manual

English.....	2
Nederlands	9
Deutsch.....	16
Français.....	23
Español.....	30
Türk.....	37
Русский.....	44

IMPORTANT: Read all instructions carefully before using this product. Retain this owner’s manual for future reference. The specifications of this product may vary from this photo and are subject to change without prior notice.

Table of contents

Safety precautions.....	2
Electrical safety.....	3
Overview drawing, parts list, and tools	3
Assembly instructions.....	5
Operating the controller	6
Workout	7
Warm up and cool down routine.....	8

Safety precautions

Read all instructions carefully before operating. Retain this owner’s manual and keep all original transaction receipts for future reference.

- Consult your physician before beginning your exercise program.
- If you experience any irregular physical conditions such as dizziness, severe muscle or joint pain or pain in the chest, stop exercising and consult your physician immediately.
- Keep children and pets away from the machine while in use.
- Do not wear loose clothing while exercising.
- Keep a minimum safety clearance during operation.
- Do not operate the exercise equipment if it is damaged.
- Set up Pedal Exerciser as directed in Assembly Instructions and place unit on a flat, stable, non-slip surface.
- Perform 3 to 5 minutes of warm-up and stretching before beginning your exercise.
- Each exercise should be performed in a controlled manner. Always start exercising slowly.
- Never exercise to exhaustion.
- Do not stand on the Pedal Exerciser.
- The Pedal Exerciser is designed as a home use exercise device.
- Users must have full control of their muscles to operate this unit. It may not be suitable for people who are quadriplegic or paraplegic. Please consult your doctor or physician.
- Exercise equipment is to be regularly inspected for maintenance for loose or broken parts.
- Have the unit serviced by a qualified service technician. Do not attempt to service the unit yourself.

Storage Precautions:

Do not store the unit in any place where it will be subject to high or low temperatures. Do not expose to direct sunlight for extended periods of time. Do not store the unit in a humid or dusty area.

Electrical safety


When using an electrical appliance, basic precautions should always be followed, including the following:

Read all instructions before using this unit.

WARNING - To reduce the risk of burns, fire, electric shock, or injury to persons:

- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts.
- Do not operate under blanket or pillow. Excessive heating may occur and cause fire, electric shock, or injury to persons.
- Close supervision is necessary when this unit is used by or near children, physically challenged or disabled persons.
- Use this unit only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this unit if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged or dropped into water. Return the appliance to a service center for examination and repair.
- Do not carry this appliance by supply cord or use cord as a handle.
- Keep the cord away from heated surface.
- Never operate the unit with the air openings blocked; keep the air openings free of lint.
- Never drop or insert any objects into any openings.
- Do not use outdoors.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- To disconnect, turn the speed to the off position, then remove plug from outlet.
- Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
- For Household Use Only.
- Power requirement: 220 V AC

Overview drawing, parts list, and tools

(24) Hexagon Socket
Pan Head Cap Bolt
6 PCS

(31) Spring Washer
6 PCS






Wrench
1 PC

Allen Wrench
1 PC

Part #	Description	Q'ty	Part #	Description	Q'ty
O01	Controller	1	O23	Phillips-Slot Pan Head Bolt M6x12	6
O02	Rear Stabilizer Ø50.8x380	1	O24	Hexagon Socket Pan Head Cap Bolt M8x12	6
O03	Front Stabilizer Ø50.8x220	1	O25	Cross Recessed Pan Head Tapping Screw ST2.9x10	2
O04	Main Frame	1	O26	Cross Recessed Pan Head Tapping Screw ST4.2x15	7
O05	Magnet Bracket Ø27.8x10	1	O27	Controller Bracket A	1
O06	Right Cover	1	O28	Washer Ø5.1xØ10.3x0.8t	1
O07	Left Cover	1	O29	Spring Washer Ø5.1xØ8.1x1.1t	1
O08	Right Crank	1	O30	Spring Washer Ø6.4xØ9.6x1.5t	6
O09	Left Crank	1	O31	Spring Washer Ø8.4xØ13.5x2.5t	6
O10	Plastic Cap	2	O32	Controller Bracket B	1
O11	PCB Support	6	O33	Right Handgrip Pedal	1
O12	Stabilizer End Cap	4	O34	Left Handgrip Pedal	1
O13	Motor	1	O35	Power Cord	1
O14	Square Sensor	1	O36	Non-slip Mat	1
O15	Square Magnet	1	O37	Power Socket	1
O16	Power Control Board	1	O38	Red Wire	1
O17	Right Foot Pedal JD-23A	1	O39	Black Wire	1
O18	Left Foot Pedal JD-23A	1	O40	EVA Cushion 30x50xT3.0	2
O19	Hexagon Head Bolt M6x40	2	O41	Rubber Bumper Pad Ø20xØ4.2x2	2
O20	Hexagon Flange Bolt M8x20	2	O42	Earth wire	1
O21	Nut M6	2	O43	Cross Recessed Pan Head Tapping Screw ST4.0x15	4
O22	Cross Recessed Pan Head Bolt M5x10	1	O44	Cross Recessed Pan Head Tapping Screw ST4.2x20	2

Assembly instructions

1. Hold the unit in an upright position. The Controller (1) side facing up and the spiral cable side is on the bottom. Then insert the Front Stabilizer (3) into the Main Frame (4).
2. Attach the Front Stabilizer (3) to the Main Frame (4) with three Hexagon Socket Pan Head Caps (24) and three Spring Washers (31). (See 'overview drawing, parts list, and tools')
3. Insert the rear stabilizer (2) into the main frame (4).
4. Attach the Rear Stabilizer (2) to the Main Frame (4) with three Hexagon Socket Pan Head Caps (24) and three Spring Washers (31). (See 'overview drawing, parts list, and tools')
5. Connect the Left Foot Pedal (18) to the Left Crank (9) assembly where would position your left foot when exercising. Thread it into the crank assembly in the counter-clockwise direction by the tool provided. (See figures 5 and 5-1). Connect the Right Foot Pedal (17) to the Right Crank (8) assembly by threading it in the clockwise direction.

Note: The left/right pedals and left/right cranks are marked with "L" and "R" (left and right)

6. Plug the indicated end of the Power Cord (35) into the socket on the machine.

Operating the controller

Button function

MODE: Press MODE button to select function display.

UP: To increase the rotational speed of the pedal. The RPM levels vary from level 0 to 5.

DOWN: To decrease the rotational speed of the pedal. The RPM levels vary from level 5 to 0.

RPM LEVEL (Rotation Speed Level)	RPM (Unloaded)
0	0 (STOP)
1	20 ~ 25
2	30 ~ 35
3	40 ~ 45
4	50 ~ 55
5	60 ~ 65

Display function

SCAN: Displays each function for 6 seconds.

TIME: Displays the time elapsed or the time remaining in your workout.

Count up: Without setting the time value, user may press the UP button on the controller to start exercise. Time starts counting up from 0:00 to 99:59 minutes per 1 second increment. The speed of pedal revolution can be adjusted by pressing the UP or DOWN button on the controller. When time counts up to 99:59, the pedals will stop to revolve automatically. The STOP sign will appear on the screen of the controller. Press UP button, you can start to exercise again.

Count down: User can set exercise time before workout. Press MODE button to select TIME function and then press and hold MODE button for 2 seconds, TIME will flash for setting exercise time. Press MODE button to set exercise time. Each increment is 1 minute. The set exercise time range is from 5:00 to 30:00 minutes. Once you set exercise time and then press UP button to start exercise, time starts counting down from the set exercise time to 0:00 per 1 second backward. The speed of pedal revolution can be adjusted by pressing the UP or DOWN button on the controller. When setting exercise time counts down to 0:00, controller will alarm for 5 seconds and then the pedals will stop to revolve automatically. The STOP sign will appear on the screen of the controller. Press UP button, you can start to exercise again.

RPM: Displays the current RPM (Revolutions per Minute) level.

COUNT: One count for each pedal revolution.

CALORIES: Estimated calories burn during exercise.

To select function display: Press MODE button.

Press and hold MODE button for 2 seconds to reset all functional values to zero when it is not in TIME function. If stop exercising for over 4 minutes, monitor will turn off.

Workout

How to use your OxyCycle 2 Pedal Exerciser

The OxyCycle 2 Pedal Exerciser can only be used after the power cord is plugged into the electrical outlet. Do NOT use the OxyCycle 2 Pedal Exerciser without the power plugged in as this will damage the equipment. When using the OxyCycle 2 Pedal Exerciser, the proper way to exercise is to set the desired rotational speed of the pedal (RPM level) by pressing the UP or DOWN button on the controller. Pedaling at excessive speed or applying excessive force will activate the overload protection circuitry to stop the motor, allowing the user to retract feet/hands. The motor will turn off for 2-3 seconds and then turn on again automatically. If the motor overload situation persists (motor turns off and then turns on) for 10 times, the motor will stop working completely. The STOP sign will flash on the controller and the controller will alarm for 8 seconds. Press the DOWN button until the RPM level reaches 0. Then you may press the UP button to start exercising again.

Leg Exercise

Place the Pedal Exerciser on a non-slip mat and sit comfortably in a chair, the device directly in front of you. The knee joint should not be flexed less than 90 degrees. Wear rubber soled shoes and adjust foot pedal straps to the proper length. Barefooted use of the OxyCycle is not recommended.

WARNING


Do not stand on the unit to exercise. Only exercise in a seated position. Make sure the seat does not move or pivot. Do not sit in a chair with casters.

NOTE: Use the foot pedals for leg exercises.

Arm Exercise

Place the unit directly in front of you on a table top (on the included anti-slip mat). Use the hand grips for arm exercises.

Warm up and cool down routine

The WARM-UP is an important part of any workout. The purpose of warming up is to prepare your body for exercise and to minimize injuries. Warm up for two to five minutes before aerobic exercising. It should begin every session to prepare your body for more strenuous exercise by heating up and stretching your muscles, increasing your circulation and pulse rate, and delivering more oxygen to your muscles.

COOL DOWN at the end of your workout, repeat these exercises to reduce soreness in tired muscles. The purpose of cooling down is to return the body to its resting state at the end of each exercise session. A proper cool-down slowly lowers your heart rate and allows blood to return to the heart.

1. Head rolls

Rotate your head to the right for one count, you should feel a stretching sensation up the left side of your neck. Then rotate your head back for one count, stretching your chin to the ceiling and letting your mouth open. Rotate your head to the left for one count, then drop your head to your chest for one count.

2. Shoulder lifts

Lift your right shoulder toward your ear for one count. Then lift your left shoulder up for one count as you lower your right shoulder.

3. Side stretches

Open your arms to the side and lift them until they are over your head. Reach your right arm as far toward the ceiling as you can for one count. Repeat this action with your left arm.

4. Quadriceps stretch

With one hand against a wall for balance, reach behind you and pull your right foot up. Bring your heel as close to your buttocks as possible. Hold for 15 counts and repeat with left foot.

5. Inner thigh stretch

Sit with the soles of your feet together and your knees pointing outward. Pull your feet as close to your groin as possible. Gently push your knees toward the floor. Hold for 15 counts.

6. Toe touches

Slowly bend forward from your waist, letting your back and shoulders relax as you stretch toward your toes. Reach as far as you can and hold for 15 counts.

7. Hamstring stretches

Extend your right leg. Rest the sole of your left foot against your right inner thigh. Stretch toward your toe as far as possible. Hold for 15 counts. Relax and then repeat with left leg.

8. Calf/achilles stretch

Lean against a wall with your right leg in front of the left and your arms forward. Keep your left leg straight and the right foot on the floor; then bend the right leg and lean forward by moving your hips toward the wall. Hold, then repeat on the other side for 15 counts.

MVS In Motion
Krommebeekpark 11
8800 Roeselare
Belgium - Europe
www.mvs-in-motion.com

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